Referencia: 032191020

Ficha Técnica

Producto: O/F Glucose Medium (Hugh y Leifson)

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Ficha Técnica

Producto: O/F Glucose Medium (Hugh y Leifson)

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- ATLAS, R.M., L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press, Inc. London.
- COWAN, S.T. (1974) Cowan and Steel's manual for the identification of medical bacteria. 2nd ed. Cambridge University Press. UK.
- DOWNES, F.P. & K. ITO (2001) Compendium of Methods for the Microbiological Examination of Foods. 4th ed. APHA. Washington. DC. USA.
- FORBES, B.A., D.F. SAHM & A.S. WEISSFELD (Eds.) (1998) Bailey & Scott's Diagnostic Microbiology. 10th ed. Mosby. St Louis. MO. USA.
- HORWITZ, W. (2000) Official Methods of Analysis of the AOAC International 17th ed. Gaithersburg. MD. USA.
- HUGH, R. & E. LEIFSON (1953) The taxonomic significance of fermentative vs. Oxidative metabolism of carbohydrates by various Gram negative bacteria. J. Bact 66:24.
- ISENBERG, H.D. (1992) Clinical Microbiology Procedures Handbook. ASM Press. Washington. DC. USA.
- ISENBERG, H.D. (1998) Essential Procedures for Clinical Microbiology. ASM Press. Washington. DC. USA.
- MURRAY, P.R., E.J. BARON, J.H. JORGENSEN, M.A. PFALLER & R.H. YOLKEN (Eds.) (2003) Manual of Clinical Microbiology. 8th ed. ASM Press. Washington. DC. USA.
- MacFADDIN, J.F. (1985) Media for Isolation-cultivation-identification-maintenance of Medical Bacteria. Vol. I. Williams & Wilkins. Baltimore. MD. USA.
- US FDA (Food and Drug Administrations) (1998) Bacteriological Analytical Manual. 8th ed. AOAC International. Gaithersburg. MD. USA.
- WINN, W., S. ALLEN, W. JANDA, E. KONEMAN, G. PROCOP, G. WOODS & P. SCHRECKENBERGER (2006) Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott, Williams & Wilkins. Philadelphia. PA. USA.
- ISO 11133:2014/ Adm 1:2018/ Adm 2:2020/ Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.