



SURFACE™



Technical Specification Sheet

SurFACE™ SPONGE-STICK

WITH BUFFERD PEPTONE WATER



Item No.	10006950
Description	Cellulose sponge attached to a quick release handle in Twirl-Tie™ Bag, individually packed & dosed with 10 ml of Buffered Peptone Water
Intended Use	For environmental surface sampling to collect and ship microbial samples
Components	Cellulose sponge, quick-release handle, 10 ml buffer, 11.5 x 30 cm ag
Buffer	BUFFERED PEPTONE WATER Neutral Peptone 10 g/L, Potassium phosphate 1.5 g/L, Sodium Chloride 5 g/L, Sodium Phosphate di-basic 3.5 g/L
Sterilization	Sterilized by irradiation at dosage of 35 KGy
Shelf Life	12 months from date of production
Storage Conditions	4-25 °C; protect from light, heat and moisture
Pack Size	100 in a box
Quality Control	Recovery of microorganisms /inhibitory effect Sterility Neutralizing effect
Conformity	Product conforms with ISO 18593:2018 "Microbiology of the food chain - Horizontal methods for surface sampling"

Date of Creation: 05. May 2022