



SteriSponge



Reference A04812

**Sponge stick, buffered peptone water
in blue sampling bag with laser pre-cut and resealable label**



**ISO 18593 standard
ISO 17604 standard**

Use:

Sterile sampling device adapted to sample carcasses and work surfaces in food industries.
Used to check the environment during production when low probability or absence of disinfectant.
Pathogens detected: Listeria, Salmonella, other ones.

Features:

- **Sponge stick**
 - Premium quality cellulose
 - Size 8 x 4 cm and 1,2 to 1.6 cm thick
 - Blue plastic made stick, with release clamps
- **Diluent**
 - Buffered peptone water (BPW)
 - Calibrated dose: 10 ml
- **Blue sampling bag**
 - Easy to open with laser pre-cut, 20 x 32,5 cm
 - Large writing area for the sample information
 - Resealable label to close the bag after sampling



Référence	Custom code	Packaging
A04812	3921140000	100 units per box

Sterility and traceability:

Sterilization	Sterilization certificate	Shelf life	Traceability*
Sterilization by gamma rays Calibrated dose minimum 20 kiloGrays	in each box	12 months	On each product/box



Before



After

* Printed information: reference, designation, batch number, diluent code, production and expiry date.

19/03/2024

