

Product Specifications

CL01.1343

IONEX

Magnesium standard solution 10000 µg/ml for IC (Plasma HIQU)

Magnesium standaard oplossing
Magnésium étalon solution
Magnesium Standard Lösung
Magnesio solución patrón
Magnesio soluzione standard

Solution contains 9900 - 10100 µg/ml of Mg

83.63 g MgCl₂·6H₂O / l H₂O

Traceability: Reference standards acc. NIST SRM 3131a

Mol weight:
Density: 1,02 g/ml
HS Nr: 28273100
CAS Nr: 7791-18-6
UN:
ADR:

Specifications: